

IN JUST THREE SECONDS ITS READY TO TRAVEL WITH YOU ANYWHERE.



Mobility with Style

ATTENDANT CONTROLLER FEATURES





Ÿ Attendant controller is a dual control device designed for attendant/caregivers to help the user while walking behind the power wheelchair.

REVO

- Ÿ The attendant controller designed for exclusively for the PWC use only.
 - Ÿ The attendant controller is attached close to hand grips.



- Ÿ The attendant controller can be connected with joystick Knob of every EVOX-WC models.
- Ÿ Provide the best possible care for attendant/caregivers while achieving normal driving speed.
- Ÿ Very easy to handle & operate.



FULLY AUTOMATIC WITH REMOTE AND JOYSTICK MODEL WC-109A

















EVOX 5LPM OXYGEN CONCENTRATOR


OXYGEN

NEBULIZE

@EVO

The first one to focus on fluctuation of an oxygen concentrator and can control the fluctuation rate within 0.1%.

Purity For Surety



Contributing to Healthy Life

Models and Functions

Model No	Max O₂ Flow	O ₂ Purity	Volts	Watts	O₂ Pressure	Nebulizer	Timer	Alarm	O ₂ 00.0% Display
EVOX-5S	5L/min	93%±3%	AC 220V	350W	≤70KPa	S	YES	YES	YES

Technical Specifications

System Map	Functions	EVOX-5S
Driving System	Oxygen Flow Oxygen Concentration Oxygen Output Pressure Oxygen Number Display	0.5⊡5L/min 93% +3% 40-70KPa YES
Nebulizer System	Nebulizer Mode Individual Nebulizer Outlet Display Method Operating Way Nebulizing Particle Size Accessories	YES YES On LED screen By touch panel button 3.6µm±25% Nebulize Kit and mask
Screen Display System	Display Material What Can Be Displayed ?	All by LED lights "Oxygen" light on "Nebulize" light on "Oxygen Percentage" light on "Each Running Time" light on "Accumulated Time" light on "Timer Setting Number"
Operating System	Operating Panel Oxygen Nebulize Timer "+" and "-"	Smart Touch Start to produce Oxygen Start to run Nebulizer Setting from 0 to 3 hours (15/30 mins increase or decrease)
Safety System	Low Oxygen Alarm Oxygen Pressure Protection Alarm Pressure Faulty Alarm Power off Alarm Compressor Faulty Alarm Overheated Protection	YES YES YES YES YES YES
Electrical System	Power Consumption Auto Power off by timer Screen Sleep Mode Noise Level	350 W YES YES <46db(A)
Packing	Machine Body Size Package Size Net Weight Gross Weight Body Material	380x270x580mm 460x350x680mm 15.6kg 19 kg ABS, Glaze

EVOX 10LPM OXYGEN CONCENTRATOR

The first one to focus on Fluctuation of an oxygen concentrator. The first one to control the Fluctuation rate within 0.1%.



Models and Functions										
Model	Max O2 Flow	O2 Purity	O2 outlet	Volts	Watts	O2 Pressure	Noise	Timer	Alarm	O2% Display
EVOX-10S	S 10L/min	93%±3%	Single	AC 220~240∨	′ 630W	≤70KPa	≤55dB(A)	YES	YES	YES
EVOX-10D) 10L/min	93%±3%	Double	AC 220~240V	630W	≤70KPa	≤55dB(A)	YES	YES	YES

Accessories							
Model	Accessories	Accessories					
Model	1 Nasal oxygen cannula						
	1 Humidifier bottle						
Evox-10S	1 Humidifier bottle intake tube	1 Humidifier bottle intake tube					
	1 User Manual						
Technical Creatification							
Technical Specificatior							
System Map	Functions	Evox-10S					
	Oxygen Flow	10 Liter Per Minute					
	Oxygen Concentration	93% ± 3%					
Oxygen Driving System	Oxygen Output Pressure	40~70 KPa					
	Oxygen Number Display	YES					
	Display Material						
		All by LED lights					
		"Oxygen" light on					
	What Can Be Displayed ?	"Purity Percentage" light on "Each Running Time" light on					
Screen Display System	What Can be Displayed ?	"Accumulated Time" light on					
		"Timer Setting Minutes"					
	Operating Panel	E- Touch					
Operating System	OXYGEN	Start to produce Oxygen					
eperating eyeteti	TIMER "+" and "-"	Setting from 0 to 3 hours					
		(15mins increase or decrease)					
	Low Oxygen Alarm	YES					
	Oxygen Pressure Protection Alarm	YES					
6 Safety System	Pressure Faulty Alarm	YES					
o curcty cystem	Power off Alarm	YES					
	Compressor Faulty Alarm	YES					
	Overheated Protection	YES					
	Machine Body Size	414*337*666 mm					
Packaging Details	Net Weight per Unit	25.7 kg					
	Operating Temperature	41°F to 113°F (5°C to 45°C)					
	Operating Humidity	30% to 80% RH					
Operating Condition	Operating Atmospheric Pressure	613-1060hp					
	Storage Temperature	14°F to 122°F(-10°C to 50°C)					
	Storage Humidity	20 to 90% RH					





Product Parameters

	Туре	Auto CPAP					
Basic Information	Device model	Noah Sleep 5A					
mormation	Working models	CPAP, Auto CPAP	PAP, Auto CPAP				
	Identification of CSA	Yes					
	Expiration pressure release	NoahSleep Elex 1-3 level					
	Auto On/Off	Yes					
	Leaking compensation	Yes					
	Manual Titration	Yes					
Functions	Therapy data review through screen	P90, AHI, SNI, LEAK90, Average pressure	CAI, AI,				
	Therapy data recording	Therapy data recording up to 5 years, reading through screen and PC software					
	PM2.5 /Pollen Filter	Folded PM2.5/Pollen filter(Unfolded area up to 20 cm2)					
	Mask leaking test	Yes					
	NoahSleep Dry Technology	30-min lower speed running of the blower keeps the device and tube dry					
	Alerts	Power down alert, leak alert, alarm clock, filter change alert, and customized alerts					
	Pressure	4-20cmH2O	Ramp	0-60min			
Main device	Dimensions	160 mm × 148 mm × 99.5 mm (Main device) 247 mm × 148mm × 99.5 mm (Humidifier included)					
parameters	weight	Main device: 0.88kg; incl. humidifier: 1.26kg					
	Noise	≤28dB	Screen	3.5 inch touch screen			
	Adapter	100-240V AC,24V DC	Language	English.			
Humidification system	Function	Heated humidification, 5 levels to adjust, infrared communication main device					
-	Humidifier volume	300ml					

BIPAP NoahSleep series

BIPAP Non-Invasive Ventilator

CHOOSE A BETTER PRODUCT FOR A HIGHER QUALITY SLEEP



@EVOX

4.3 inch



4.3 Inches high resolution colorful screen.

Exquisite UI design, multi-theme of the monitored parameters on one screen.

Touchscreen and high flexible navigator knob.

Easy to operate.

Light-Weighted 1.26kg. Convenient for home and traveling.

NoahSiceptes

Silent

The noise level achieves less than 30dB contributed by the research from experienced team.

Portable


Detachable design of the main unit and humidifier meets different demands of the user.

ACCESSORIES


CPAP: Continuous Positive Airway Pressure mode, suitable for obstructive Sleep Apnea.

- AUTOS: Intelligent regulation of inspiratory pressure to ensure a more smooth and comfortable breath.
 - S: Spontaneous ventilation mode, suitable for patient with spontaneous respiratory
 - T: Timed mode, suitable for patients without spontaneous respiratory capacity.
 - **S/T:** The machine automatically switches between S & T modes according to patient's respiratory condition.
 - VGPS: Volume Guarantee Pressure Support, meets patient's tidal volume requirement with its minimum pressure.

PRODUCT SPECIFICATION							
Product Model	NoahSleep25A	NoahSleep25ST	NoahSleep30ST				
Ventilation Modes	CPAP, S, Auto S	CPAP, S, T, ST, VGPS	CPAP, S, T, ST, VGPS				
ΙΡΑΡ	6-25hPa	6-25hPa	6-30hPa				
EPAP	4-23hPa	4-23hPa	4-25hPa				
MIN IPAP	6.0-25.0hPa	6.0-25.0hPa	6-30hPa				
ΜΑΧΙΡΑΡ	6.0-25.0hPa	6.0-25.0hPa	6-30hPa				
Target Tidal Volume	/	200ml-1500ml	200ml-1500ml				
Res. Rate	10BPM	10BPM / 5-40BPM	10BPM / 5-40BPM				
Ti (sec.)	2s	2s	2s				
	(Backup frequency)	(Backup frequency) $/ 0.3$ -4s	(Backup frequency) / 0.3-4s				
Ti min (sec.)	0.5-4s	0.5-4s	0.5-4s				
Ti max (sec.)	1-4s	1-4s	1-4s				
Trigger and Cycle	1-6	1-6	1-6				
Rise Slope	1-6 1-6		1-6				
Monitoring Parameters	Pressure, Tidal volume, Minute ventilation, Respiratory rate, Leakage						
Oximetry Accessory	Optional						
Therapy Data Storage	Yes, SD card						
NoahSleep -Dry Drying Technology	Yes						
Auto On/Off	Yes						
Noise Level	<30 dB						
Screen Size	4.3 Inches Color Touchscreen						
Weight	Main unit: 0.82kg With humidifier: 1.26kg						
Adaptor	AC input: 100-240V, DC output: 24V						



Pressure rise slope: 6 levels adjustable



VGPS (Volume Guarantee Pressure Support) mode, ensures the ventilation volume target.

TWINPAP series COMBINATION OF AUTOCPAP & AUTOBIPAP



CPAP/BIPAP

For Sleep Apnea Hypopnea Syndrome(SAHS) NIV FULL-FACE MASK

田 Technical Data 日本				
	CPAP	TWINPAP B-19		
	C-09			
Modes				
Μ	СРАР	СРАР		
0	Auto CPAP	Auto CPAP		
D		S		
E		Т		
S		S/T		
		Auto BIPAP		
Parameters				
Pressure Range	4-20 cmH2O	4-25cmH2O		
Belex	2-4cmH2O	2-4cmH2O		
Ramp	0-60 min	0-60 min		
Esens (1-6 levels)		\checkmark		
Isens (1-6levels)				
Auto-leak Compensation	\checkmark			
Auto-altitude Adjustment	\checkmark			
Auto-On/Off	\checkmark	\checkmark		
Pressure Titration	\checkmark			
Noise Level (10cmH2O)	<30dBA	<35dBA		
Heating Humidifier				
Adjustable Humidify	\checkmark			
Replaceable Water Chamber	\checkmark			
Data Download				
Data Details	\checkmark			
Data Summary	\checkmark			
Download to PC	\checkmark	\checkmark		

NoahSleep series





- Plot No. EL-635, Phase-9 Industrial Area Sector-66A, SAS Nagar (Mohali) Punjab-160062
- 🔇 오 7696052635, 0172-5099774
- 🔀 info@evoxcare.in
- www.evoxworld.com

